CE  
0086

**Instructions for Use**  
**ClearView MV™ Atrial Depressor**  
**Catalog No. 400-005**

**STERILE**

## **READ ALL INSTRUCTIONS BEFORE USE!**

**Federal (USA) law restricts this device to sale by or on the order of a physician.  
Single Use Device. Do not resterilize.**

### **PRODUCT DESCRIPTION**

The ESTECH ClearView MV™ Atrial Depressor is a single-use sterile device. The atrial depressor provides exposure of the mitral valve via a plastic paddle and malleable shaft for retraction of the posterior atrial wall. It has an integrated suction channel for clearing the atrium of fluid during cardiac surgery. The atrial depressor is designed for use with the ESTECH SI Retractor System (ref 400-400).

### **INDICATION FOR USE**


The ESTECH ClearView MV™ Atrial Depressor is intended for use in cardiac surgical procedures to aid in exposure of the mitral valve.

### **WARNINGS AND PRECAUTIONS**

- Federal Law (USA) restricts this device to sale by or on the order of a physician.
- Do not use this device if package is damaged or opened.
- Do not use this device if it shows signs of damage, e.g., crimps, kinks, or crushed areas.
- This device is for single use only. **DO NOT REUSE. DO NOT RESTERILIZE.**
- Proper surgical procedures and techniques are necessarily the responsibility of the medical professional. These instructions are furnished for informational purposes only. Each surgeon must evaluate the appropriate use of this device, case by case, based on medical training, experience, and the type of surgical procedure employed.
- Dispose of used products in accordance with established hospital protocol for hazardous waste.
- Verify all user connections of this device are secure prior to use.

### **USAGE INSTRUCTIONS**

1. Using aseptic technique, remove the atrial depressor from sterile packaging and place in the sterile field.
2. After the SI Retractor (ref pn 400-400) has been installed and the atrium has been opened, place the paddle end of the depressor into the atrium with the suction inlet facing down.
3. Loosely attach the mount screw at the proximal end of the depressor to the threaded hole at the left side of the SI Retractor.
4. Position the depressor as desired and tighten down the mount screw to hold the atrial depressor firmly in place.
5. When the mount screw is tight, manipulate the malleable shaft to the desired final position for depressing the atrial wall.
6. Connect the proximal slip connector to the return line of the cardio pulmonary bypass circuit.
7. When the procedure is finished, disconnect the proximal connector.
8. Remove the mount screw and discard the entire atrial depressor assembly.

**Figure 1**


**WARRANTY**

AtriCure provides no warranties, either express or implied, with respect to the product, including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. AtriCure shall not be liable to any person or entity for medical expenses or any direct, incidental or consequential damages caused by any use, defect, failure or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort or otherwise. No person has any authority to bind AtriCure to any representation or warranty with respect to the product.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this DISCLAIMER OF WARRANTY is held to be illegal, unenforceable or in conflict with applicable law, by a court of competent jurisdiction, the validity of the remaining portions of this DISCLAIMER OF WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this DISCLAIMER OF WARRANTY did not contain the particular part or term held to be invalid.

Manufactured by:  
**AtriCure Incorporated**  
7555 Innovation Way  
Mason, Ohio 45040 USA  
Customer Service:  
1-866-349-2342 (toll free)  
1-513-755-4100 (phone)

European Representative:  
Köntges SPRL  
Avenue Helleveldt 35  
B-1180 Brussels Belgium  
Tel: +32 (0) 2 375 51 63  
FAX: +32 (0) 2 375 89 06  
e-mail: herbert.kontges@skynet.be

CE  
0086

**Graphic Symbols for Device Labeling**

<p>Catalogue Number</p>	<p>Batch Code</p>	<p>Quantity</p>	<p>Date of Manufacture</p>
<p>Use by date</p>	<p>Manufacturer/Company Address</p>	<p>Authorized Representative in the European Union</p>	<p>Sterilized by ethylene oxide treatment</p>
<p>Does not contain Natural Rubber Latex</p>	<p>Caution: U.S. Federal law restricts this device to sale by or on the order of a physician or other licensed practitioner.</p> <p>Rx ONLY</p>	<p>Consult Instructions For Use</p>	<p>Do not reuse</p>
<p>CE mark, Product conforms with the essential requirements in the European Medical Devices Directive 93/42/EEC.</p>	<p>Storage temperature range</p>	<p>Storage humidity range</p>	

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