Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

COBRA AFfirm[®] Bipolar Pacing Probe

Catalog Number 600-003

U.S. Patent No.'s 7,542,807; 7,957,820; 8,096,990; Other patents pending.

Instructions for Use

TABLE OF CONTENTS


ENGLISH.....	2
NEDERLANDS.....	6
FRANÇAIS.....	11
DEUTSCH.....	16
ITALIANO.....	23
ESPAÑOL.....	28
HRVATSKI.....	33
SLOVENČINA.....	38

Sterile – Single Use Only

AtriCure, Incorporated
6217 Centre Park Drive
West Chester, Ohio 45069-3866
USA
Customer Service Tel: 866.349.2342
Email: customerservice@AtriCure.com
www.AtriCure.com

Köntges SPRL
Avenue Hellevelt 35
B-1180 Brussels Belgium
Tel: +32 (0) 2 375 51 63
Fax: +32 (0) 2 375 89 06
Email: herbert.kontges@skynet.be

ENGLISH

Instructions for Use

Carefully read all instructions prior to use. Observe all contraindications, precautions, cautions and warnings noted in these directions. Failure to do so may result in patient complications.

AtriCure relies on the physician to determine, assess and communicate to each patient all foreseeable risks of the procedure.

Product Description

The AFfirm Bipolar Pacing Probe is a pacing and sensing device which provides the ergonomic convenience of a handle and malleable shaft combined with an atraumatic bipolar electrode tip that is used epicardially and does not require temporary implantation in heart tissue.

Intended Use

The AFfirm Bipolar Pacing Probe is intended to be used within the indications for use for the Cobra Surgical System. The Electrode is designed to provide transient cardiac pacing or recording for the assessment of electrical isolation / conduction block of ablation lesions in the surgical treatment of arrhythmias.

Device Indication

Temporary cardiac stimulation or ECG monitoring.

Inspection Prior to Use

- Carefully inspect the package prior to use for any breach of the sterile barrier or damage to the contents.
- Inspect device for physical damage, including electrical insulation on the cables and the catheter shaft.
- If the sterile barrier integrity is compromised or the contents are damaged, DO NOT USE and contact your AtriCure representative. Note: Unopened, undamaged package is guaranteed sterile up to "Use Before" date.

Sterilization Method

The sterilization process for the device is gamma radiation and provides sterility assurance level 10^{-6} in compliance with the obligatory requirements of ISO 11137.

Precautions

- The use of this device should be performed by an experienced physician only.
- The device is for single use only. Do not reuse, reprocess or re-sterilize. Reprocessing or re-sterilization may compromise the structural integrity of the device and/or lead to device failure which in turn may result in patient injury, illness or death. Reuse, reprocessing or re-sterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. This device is for one time use only!
- Do not use if sterile barrier is damaged.
- Replace damaged equipment. If damage is found call your AtriCure representative.

- Defibrillation and life support equipment and supplies should be placed on standby and be kept immediately available during device manipulation, pulse generator connection, and removal.

Cautions

- Bare electrode tips and lead wire connectors constitute a direct path to the heart. Small currents (approximately 10 μ amps) may be sufficient to cause fibrillation. Only trained personnel should handle the AFfirm Bipolar Pacing Probe and connecting pins.

Warnings

- The device(s) should be used by physicians thoroughly trained in the techniques of invasive surgical procedures and in the specific approach to be used.
- Line powered equipment used in the vicinity of the patient must be properly grounded.
- Avoid the possibility of unintentional contact between the device tip electrodes or pacing lead connectors and any conductive surface contact. The pacing lead connectors should only be connected to a temporary pacemaker or to a recorder designed for safely recording biopotentials.

Contraindications

DO NOT USE WHEN A PACEMAKER HAS BEEN IMPLANTED.

Potential Complications

The following potential risks or discomforts may be associated with cardiac surgical procedures. The frequency and severity of these events can vary, and may necessitate additional medical intervention, including additional surgery. Strict adherence to the following instructions before use will help reduce the incidence of complications.

Allergic reaction, Arrhythmias, Cardiac or respiratory arrest, Cardiac valve damage, Chest pain, Death, Embolus, air embolus, Hematoma / ecchymosis, Hemorrhage, Infection, Perforation, Pericardial effusion, Pericarditis / pleuritis, Pseudoaneurysm, Pulmonary edema, Stroke, Tamponade, Vasovagal reaction.

DIRECTIONS FOR USE

Pacing:

1. Use a compatible extension cable to connect the two proximal connector ends of the AFfirm Bipolar Pacing Probe to a temporary pacemaker.
2. Set the Output Pulse Amplitude (V) to the desired level on the temporary pacemaker.
3. Set the Output Pulse Duration (ms) to the desired level on the temporary pacemaker.
4. Set the Output Pulse Frequency (BPM) to the desired level on the temporary pacemaker.
5. Contact the distal bipolar electrode tip of the AFfirm Bipolar Pacing Probe to the desired target area of the epicardium (atrial or ventricular) for pacing. It is important that both electrodes contact the epicardium prior to pacing.
6. Activate the temporary pacemaker to compare physiological response to pacing stimuli for various desired target anatomical areas of the myocardium (atrial or ventricular).

Sensing:

1. Use a compatible extension cable to connect one of the proximal connectors to an ECG Recorder.
2. Contact the distal bipolar electrode tip of the AFfirm Bipolar Pacing Probe to the desired target area of the epicardium (atrial or ventricular) for ECG sensing.
3. Electrogram measurements should be done only with recorders having inputs isolated from earth ground. Non-conformance will increase the danger of fibrillation.

Note: The physician may elect to vary the procedure in accordance with judgment and experience.

Handling and storage

The package and the product tolerate normal handling, but excessive mechanical stress could cause damage to both packaging and the product. The package should be stored in a dry place at temperatures -10°C (14°F) and 60°C (140°F). Always remove the product from the inner package in a sterile environment. Before opening the package, ensure no accidental rupture has occurred. Take extreme care while removing the product. Do not use sharp instruments in handling the product to prevent insulation damage.

How supplied

The AFfirm Bipolar Pacing Probe is a single use disposable instrument which is provided sterile in individual packaging.

Contents

The packaging configuration consists of three (3) or (1) per box.

One (1) Instructions for Use per box.

End of life

Refer to the package label for the product shelf life. The shelf life indicated on the package label is valid as long as the sterility barrier is not compromised in any way. If the product's performance fails while still under warranty, please contact your local AtriCure sales representative or customer service for a RMA number to return the unit.

Disposal


After use, dispose of product and packaging in accordance with hospital, administrative and/or local, state, federal and international laws and regulations.

For technical support contact:

AtriCure, Incorporated

6217 Centre Park Drive
West Chester, Ohio 45069-3866
USA
Customer Service Tel: +1.866.349.2342
Email: customerservice@AtriCure.com
www.AtriCure.com

Köntges SPRL
Avenue Hellevelt 35
B-1180 Brussels Belgium
Tel: +32 (0) 2 375 51 63
Fax: +32 (0) 2 375 89 06
Email: herbert.kontges@skynet.be

Warranty and Limitations

AtriCure warrants that reasonable care has been used in the design and manufacture of this instrument. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular use. Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond AtriCure's control directly affect the instrument and the result obtained from its use. AtriCure's obligation under this warranty is limited to the repair or replacement of this instrument and AtriCure shall not be liable for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of this instrument. AtriCure neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. AtriCure assumes no liability with respect to instruments reused, reprocessed or re-sterilized and makes no warranties, expressed or implied, including but not limited to merchantability or fitness for intended use, with respect to such instrument.

Graphic Symbols for this Device Labeling

Catalogue Number 	Caution: See Instructions for Use 	Manufacturer/Company Address 	Sterilized using Gamma Radiation 
Batch Code 	Single Use Only 	Expiration Date 	Authorized Representative in the European Union 
Quantity 	Date of Manufacture 	Non-Pyrogenic 	Does not contain Natural Rubber Latex 
Contains no di (2-ethylhexyl) phthalate (DEHP) released from polyvinyl chloride (PVC) 	Consult Instructions For Use 	CE Mark with Notified Body - class II and III medical devices comply with EC Directive 93/42/EEC 	Caution: U.S. Federal law restricts this device to sale by or on the order of a physician or other licensed practitioner." 

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