**AtriCure® Isolator® Surgical Ablation
System
Instructions for Use**

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AtriCure® Isolator® Surgical Ablation System Instructions for Use

DESCRIPTION

The ATRICURE ISOLATOR Surgical Ablation System is comprised of the Ablation and Sensing Unit (ASU), an ATRICURE ISOLATOR device, and a footswitch. The ISOLATOR is a single patient use electrosurgical instrument designed for use only with the ASU. The ISOLATOR is used for soft tissue ablation. When activated, the ASU delivers radiofrequency (RF) energy to the linear electrodes on the insulated jaws of the ISOLATOR. The Operator controls the application of this RF energy by pressing the Footswitch.

All ISOLATOR devices are configured as vascular clamps and feature clamping jaws of various lengths and curvatures.

There are three families of ISOLATOR clamps available for use with the system.

The ATRICURE ISOLATOR (See Figure 1) Clamps feature an in-line handle with syringe-type actuation and button release mechanisms. Some ISOLATOR Devices have an integral attachment tip and related Glidepath™ Tape Instrument Guide (Guide) as device options (See Figure 1 and Figure 2). The Guide is packaged with ISOLATOR devices that have the Attachment Tip.

The ISOLATOR Synergy™ (See Figure 1) Clamps feature two pairs of opposing dual electrodes, an in-line handle with syringe-type actuation and button release mechanisms. The ISOLATOR Synergy requires the use of the AtriCure Switch Matrix (ASB3, US and EU, respectively)

There are two families of GLIDEPATH Tape Instrument Guides for use with ATRICURE ISOLATOR Devices. The GLIDEPATH Tape Instrument Guide is a single patient, surgical device designed to facilitate the guidance of surgical instruments through soft tissue during general surgical procedures


The GLIDEPATH Tape Instrument *Snap* Guide (See Figure 2) is designed to fit devices that have the integral attachment tip.

The GLIDEPATH Tape Instrument *Press* Guide (See Figure 2) is designed to fit instruments without integral attachment tips.

NOTE: Please refer to the ATRICURE ASU and ASB3 Instructions for Use for information specific to the ASU and ASB3.

ATRICURE ISOLATOR ILLUSTRATION AND NOMENCLATURE

(Figure 1)


(ATRICURE ISOLATOR INLINE SYRINGE GRIP CLAMPS)

- | | |
|--------------------|-------------------------------|
| 1. Distal Jaw | 6. Handle |
| 2. Proximal Jaw | 7. Closure Lever |
| 3. Attachment Tip* | 8. Release Mechanism |
| 4. Electrode | 9. Connector |
| 5. Jaw Heel | 10. Connector Alignment Arrow |

**Only on devices packaged with GLIDEPATH Tape Instrument Guide*

GLIDEPATH TAPE INSTRUMENT GUIDE ILLUSTRATION AND NOMENCLATURE

(Figure 2)

(Guide)

1. Lateral Tab
2. Distal Pocket
3. Instrument Attachment Tip
4. Red Elastic Leader
5. Clear Ribbon
6. Accessory Ribbon

INDICATION FOR USE

The ATRICURE Bipolar (Transpolar) System is intended to ablate soft tissue during General surgical procedures.

CONTRAINDICATIONS

The Bipolar (Transpolar) System is not indicated for contraceptive coagulation of the fallopian tubes.

Potential Complications

Device

Possible complications related to the creation of the linear lesions in cardiac tissue using a clamp-type device maybe be included but not limited to:

- Tissue Cutting
- Perioperative heart rhythm disturbance (atrial and/or ventricular)
- Postoperative embolic complications
- Pericardial effusion or tamponade
- Injury to the great vessels
- Valve leaflet damage
- Conduction disturbances (SA/AV node)
- Acute ischemic myocardial event
- Injury to unintended surrounding tissue structures, including tears and punctures.
- Bleeding requiring intervention to repair.
- Extension of cardiopulmonary bypass

Procedure

Serious adverse events that may be associated with surgical ablation procedures on the heart (stand alone or concomitant to other cardiac surgery), include:

- Death
- Excessive bleeding related to the procedure (defined as bleeding which requires >3 units of blood products and/or surgical intervention),
- Cardiac tamponade (if either open or catheter drainage is required),
- Pulmonary vein stenosis,
- Restrictive (constrictive) pericarditis,

- Endocarditis,
- Myocardial infarction (MI) per ACC guidelines,
- Stroke (resulting in permanent neurological deficit),
- Transient Ischemic Attack (TIA),
- Thromboembolism
- Diaphragmatic paralysis,
- Esophageal-LA fistula or esophageal rupture,
- Atrial perforation or rupture,
- Ventricular perforation or rupture,
- Atelectasis,
- Pneumonia,
- Congestive Heart Failure,
- Cardiac Valve Injury,
- Persistent Pneumothorax (requiring intervention),
- Excessive Pain and Discomfort,
- Deep Sternal Wound Infection,
- Ventricular Arrhythmia (V. Tachycardia or V. Fibrillation)
- New Sinus Node Dysfunction, and
- Drug Reaction

WARNINGS

- Do not touch the electrodes of the ISOLATOR while activating the ASU. Touching the ISOLATOR electrodes during ASU activation could result in an electrical shock or burn to the operator.
- Do not touch the electrodes of the ISOLATOR to metal staples or clips, or to sutures while activating the ASU. This may damage the ISOLATOR or tissue, or result in an incomplete ablation.
- Do not use abrasive cleaners or electro-surgical tip cleaners to clean debris from the Jaws. Use of abrasive cleaners or electro-surgical tip cleaners can damage the electrodes and result in device failure. Use saline-soaked gauze to clean debris off the electrodes.
- Do not immerse any part of the ISOLATOR in liquids as this may damage the device.

- Always wear the appropriate surgical gloves when using the ATRICURE ISOLATOR Surgical Ablation System to avoid shock/burn hazards.
- Inspect the product packaging prior to opening to ensure that the sterility barrier is not breached. If the sterility barrier is breached, do not use the ISOLATOR to avoid the risk of patient infection.
- Electrosurgery should be used with caution in the presence of internal or external pacemakers. Interference produced with the use of electrosurgical devices can cause devices such as a pacemaker to enter an asynchronous mode or can block the pacemaker entirely. Consult the pacemaker manufacturer or hospital Cardiology department for further information when use of electrosurgical appliances is planned in patients with cardiac pacemakers.
- Use ISOLATOR only for soft tissue ablation. Variations in specific procedures may occur due to individual physician techniques and patient anatomy.
- Do not drop or toss the ISOLATOR as this may damage the device. If the ISOLATOR is dropped, do not use. Replace with a new ISOLATOR.
- Do not use the ISOLATOR in the presence of flammable materials.
- Do not re-sterilize or reuse the ISOLATOR. Reuse can cause patient injury and/or the communication of infectious disease(s) from one patient to another.
- Keep the Jaws of the ISOLATOR clean of debris during surgery to avoid loss of power.

PRECAUTIONS

- Read all instructions carefully for the ATRICURE ISOLATOR Surgical Ablation System, prior to using the device. Failure to properly follow instructions may lead to electrical or thermal injury and may result in improper functioning of the device.
- Use of the ISOLATOR should be limited to properly trained and qualified medical personnel.
- Do not use of the ISOLATOR with another manufacturer's generator to avoid damage to the device, which may result in patient injury. The ISOLATOR is only compatible with the ATRICURE ASU and ASB3.
- Do not ablate tissue greater than 10 mm thick with the ISOLATOR. Tissues greater than 10 mm thick may not be fully ablated.
- Do not use the ISOLATOR for coagulation or ablation of veins or arteries.
- Inspect the area between the Jaws of the ISOLATOR for foreign matter before activating the ASU or ASB3. Foreign matter captured between the Jaws will adversely affect the ablation.

- Do not insert excessive tissue into the Jaw heel as it may result in poor ablation at the Jaw Heel.
- Do not ablate in pool of blood or other fluids as this may extend the ablation time. Users should suction excess fluids away from the jaws prior to ablation.
- Do not attempt to use an ISOLATOR that has reached its time limit expiration. The ISOLATOR has an 8-hour useful life that is tracked by the ASU. The ISOLATOR will no longer function after 8 hours of use and the ASU will display a message indicating that the ISOLATOR must be replaced.
- Do not use the ISOLATOR if signs of damaged wire insulation are noted upon inspection of the area around the Jaw heel as it may adversely affect ablation performance.
- When the ASU (RF generator) and Handpiece are used on a patient simultaneously with physiological monitoring equipment, ensure that the monitoring electrodes are placed as far as possible from the surgical electrodes. Be sure to position the Handpiece cables so that they do not come in contact with the patient or the other leads.
- Needle monitoring electrodes are not recommended for use when operating the ASU (RF generator) and Handpiece.
- Monitoring systems that incorporate high frequency current-limiting devices are recommended for use with the ASU (RF generator) and Handpiece.
- When the ASU (RF generator) is activated in conjunction with the Handpiece, the conducted and radiated electrical fields may interfere with other electrical medical equipment. Refer to the ASU IFU for more information regarding potential electromagnetic or other interference, and advice regarding avoidance of such interference.

INSTRUCTIONS FOR USE

SET UP

1. Examine the packaging of the devices to ensure the sterility of the product has not been breached. Remove the sterilized instruments from their package per standard sterile technique.
2. With the Connector Alignment Arrow symbol in the 12 o'clock position, push the Connector into the appropriate ISOLATOR receptacle on the front of the ASU or ASB3. Each ISOLATOR has a unique receptacle on the ASB3. To ensure device performance, verify proper connections to the ASB3 by consulting the ASB3 package insert. Verify that the connections between the ISOLATOR and the ASU or ASB3 are secure. If the connections are loose, do not use the ISOLATOR. Inspect the Cable and do not use the ISOLATOR if the cable is frayed or the insulation is damaged.

Insert with Connector Alignment Arrow at the 12 o'clock position


3. If the ISOLATOR is to be used with a supplied GLIDEPATH Instrument Guide (Integrated Guide – Figure 2), go to step 4.
If the ISOLATOR is not to be used with the instrument guide, go to step 20.

POSITIONING OF GLIDEPATH TAPE GUIDE

4. The Guide may be used with commercially available general dissection or surgical clamping tools (auxiliary tools) to create and maintain a dissection plane that facilitates placement of the ISOLATOR clamp

around structures per standard surgical technique.


5. Examine the Guide package to ensure the sterility of the product has not been breached. Remove the Guide from its packaging per standard sterile technique.
6. Secure the proximal end of the Guide to the sterile drape near the surgical site.
7. Insert the distal end of the auxiliary tool completely into the distal pocket of the clear ribbon portion of the Guide.

8. Maintain attachment of the distal portion of the Guide to the auxiliary tool during positioning of the Guide.

Note: Lubrication may be applied to the Guide at the user's discretion.

9. Once the desired placement of the Guide is achieved, use a grasping device to grasp one of the Lateral Tabs (Figure 2) on the Guide and remove the Guide from the auxiliary tool. Externally secure the distal end of the Guide near the surgical site.


Note: If an articulating dissection tool is used, un-articulating the device may facilitate removal of the Guide.

10. If desired, the Guide can be used for soft tissue retraction or to introduce additional Instruments through the previously created positioning plane.
11. If the Guide incorporates a snap feature, refer to steps 13-14 for instrument exchange.
12. If the Guide incorporates a press feature, refer to steps 15-19 for instrument exchange.

INSTRUMENT EXCHANGE (SNAP GUIDE)

13. Prior to attaching the Guide to the ISOLATOR, unclamp the proximal end of the leader from the sterile drape.
14. While holding the Guide as shown, in the illustration below, insert the Instrument Attachment Pin (Figure 2) into the ISOLATOR attachment tip (Figure 1). Once Guide is attached, do not attempt to remove by forcibly pulling on Guide.

INSTRUMENT EXCHANGE (PRESS GUIDE)

15. If using an AtriCure Instrument Press Guide, attach the guide to the distal tip of the ISOLATOR per standard surgical technique.
16. Use the guide to facilitate the placement of the ISOLATOR in the previously created positioning plane.
17. Carefully remove the guide from the distal jaw after ISOLATOR placement.

NOTE: The Press Guide is to be removed prior to ablation. (Refer to Step 28)

NOTE: The accessory ribbon allows the surgeon to create a dissection plane in one direction (inferior or superior) and maneuver the tape and leader through the dissection plane so that an instrument (clamp) can be used from the opposite direction.

18. If it is required to reverse the direction of device placement, the accessory ribbon is attached to the distal pocket of the primary transfer tape.
19. After creation of the dissection plane, pull the distal end of the primary tape while providing counter traction on the proximal end of the accessory ribbon so that the leader is pulled through the dissection plane.

ABLATION

NOTE: A minimum tissue incision of 12 mm is recommended for insertion of the ISOLATOR.

20. Place the targeted tissue between the Distal and Proximal Jaws.
21. Squeeze the Closure Lever to close the Jaws. Ensure that no target tissue extends beyond the Indicator Line on either the Distal or Proximal Jaws or into the Jaw Heel.
22. Activate the ASU by depressing the footswitch.

When the ASU is activated, the ASU will emit an audible tone indicating that current is flowing between the Jaws of the ISOLATOR.

When the continuous tone switches to intermittent, release the footswitch.

23. The ATRICURE ISOLATOR Surgical Ablation System measures tissue impedance and temperature throughout the ablation cycle and uses this information to control the application of energy to the tissue. The amount of energy delivered to the tissue is driven solely by tissue impedance. The System determines the minimum energy delivery required to create a transmural (full thickness) lesion based on tissue impedance and delivers only that amount of energy to the tissue. Energy delivery changes throughout the ablation cycle as tissue impedance changes. The lesion is visible as a white coloration of the tissue. The device is designed such that the lesions will not spread beyond the jaw width.

Note: Some ISOLATOR devices include a temperature probe that measures the surface temperature of the tissue being ablated 1.3 mm laterally from the electrode within the insulated jaw. An audible tone sounds when the tissue temperature reaches 70°C and the system automatically shuts off energy delivery when the tissue temperature reaches 75°C. All of the clamps, with and without a temperature probe, have been designed to maintain less than 50°C temperature outside of the clamped region.

Note: The time necessary to create a transmural lesion depends on tissue thickness, composition, and the length of tissue captured between the electrodes. The following table describes the average expected time (seconds) and energy delivery (joules) for respective tissue thicknesses. Values are expressed per unit volume of tissue captured between the electrodes. These data were obtained during ablations on ex vivo (excised bovine) tissues and will generally be lower on in vivo (live human) tissues.

Table1. Average Time vs. Energy Delivery

| Tissue Thickness | Time to Transmurality per unit volume (sec/mm ³) | | Energy Delivered per Unit Volume (J/mm ³) * | |
|------------------|--|--------|---|--------|
| Isolator | | | | |
| | AVG | STDDEV | AVG | STDDEV |
| 2 mm | 0.041 | 0.014 | 0.53 | 0.16 |
| 5 mm | 0.034 | 0.013 | 0.43 | 0.15 |
| 10 mm | 0.040 | 0.011 | 0.50 | 0.14 |
| Synergy | | | | |
| 2 mm | 0.049 | 0.007 | 0.76 | 0.11 |
| 5 mm | 0.033 | 0.006 | 0.57 | 0.10 |
| 10 mm | 0.032 | 0.009 | 0.55 | 0.16 |

*Energy Delivery per unit volume of tissue ablated is comparable between the Isolator and Synergy Product Lines and below the threshold of 0.94 J/mm³ for 2 mm tissue thickness reported for other similar commercially available ablation devices.

24. To open the Jaws, press the Release Mechanism and slowly release the Closure Lever. Do not allow the Jaws to spring back. Be aware of any surrounding tissues that could be damaged as the Jaws open.
25. Inspect the surgical area to ensure adequate ablation.
26. Between ablations, wipe the Jaws clean with a saline-soaked gauze pad. Important: For optimal performance, keep the ISOLATOR electrodes clear of coagulum. To ensure the electrodes are clear of coagulum:
 - Use a saline soaked gauze to clean the electrodes after each ablation. If coagulum is present, it is much easier to remove within the first several seconds after ablation. In a brief period of time, the coagulum could dry out making removal more difficult.
 - Check both electrodes before each ablation to ensure that the gold of the electrode is visible and coagulum is removed.
 - If the ISOLATOR is idle between ablations, clamp the jaws onto saline soaked gauze to prevent any coagulum on the electrodes from drying.
27. Repeat the ablation process as necessary.

REMOVAL AND DISPOSAL

28. To remove the Press Guide Attachment Tip from the Instrument, place a grasping instrument at the distal end of the Instrument and carefully back the Attachment Tip off the Instrument Jaw using a rotating motion.
29. Ensure Guide is removed from the surgical field prior to the completion of the surgical procedure. Discard the Guide after use. Follow local governing ordinances and recycling plans regarding disposal or recycling of device components.
30. Discard the ISOLATOR after use. Follow local governing ordinances and recycling plans regarding disposal or recycling of device components.

HOW SUPPLIED

The ISOLATOR and Instrument Guide are supplied as STERILE instruments and are for single patient use only.

Sterility is guaranteed unless the package is opened or damaged. Do not resterilize.

RETURN OF USED PRODUCT

If for any reason these products must be returned to ATRICURE, a return goods authorization (RGA) number is required from ATRICURE prior to shipping.

If the products have been in contact with blood or body fluids, they must be thoroughly cleaned and disinfected before packing. They should be shipped in either the original carton or an equivalent carton, to prevent damage during shipment; and they should be properly labeled with an RGA number and an indication of the biohazardous nature of the contents of shipment.

Instructions for cleaning and materials, including appropriate shipping containers, proper labeling, and an RGA number may be obtained from ATRICURE, Inc.

CAUTION: It is the responsibility of the health care institution to adequately prepare and identify the products for shipment.

DISCLAIMER STATEMENTS

Users must assume responsibility for approving the condition of these products before they are used. ATRICURE, Inc. cannot be held liable for any consequential damage, personal injury or damage to property nor for the misuse of these products.

If the products have been in contact with blood or body fluids, they must be thoroughly cleaned and disinfected before packing. The products should be shipped in either the original carton or an equivalent carton, to prevent damage during shipment; and they should be properly labeled with an RGA number and an indication of the biohazardous nature of the contents of shipment.

ATRICURE, Inc. will not be liable for any damage caused by the reuse of these products.

This Instruction for Use describes the procedures for proper use of the products. Any deviation from these procedures, which may compromise the function of the products, is the responsibility of the user.

| | |
|--|---|
|  <p>Non-Pyrogenic</p> | <p>Rx ONLY</p> <p>Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician</p> |
|  <p>Sterilized by Ethylene Oxide</p> |  <p>Follow instructions for use</p> |
|  <p>Single Use Only</p> |  <p>CE Mark with Notified Body - class II and III medical devices comply with EC Directive 93/42/EEC</p> |
|  <p>Expiration Date</p> |  <p>Does not contain Natural Rubber Latex</p> |
|  <p>Lot Number</p> |  <p>Manufacturer</p> |

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