Instructions for Reprocessing of Reusable Devices


TABLE OF CONTENTS

English	2
Dutch	6
French.....	10
German	14
Italian	19
Spanish.....	23

Manufactured by:
AtriCure Incorporated
7555 Innovation Way
Mason, Ohio 45040 USA
Customer Service:
1-866-349-2342 (toll free)
1-513-755-4100 (phone)

European Representative:
Köntges SPRL
Avenue Hellevelt 35
B-1180 Brussels Belgium
Tel: +32 (0) 2 375 51 63
FAX: +32 (0) 2 375 89 06
e-mail: herbert.kontges@skynet.be

ENGLISH

INTRODUCTION

The following instructions are applicable for reusable medical devices manufactured and supplied by ESTECH, unless stated otherwise in the IFU or other applicable documents. This instruction is intended for personnel with the required knowledge and training for the reprocessing of reusable devices to ensure instruments are properly cleaned and viewed as hygienically safe for use.

ESTECH recommends using the following web-sites for additional information on the reprocessing of surgical equipment.
www.a-k-i.org www.cdc.gov www.rki.de

GENERAL WARNINGS

1. During the reprocessing, always handle ESTECH instruments with care. It is mandatory to wear personal protective equipment i.e., clothing, gloves, eyewear, etc., in accordance with applicable national and local health & safety procedures.
2. ESTECH does not recommend any disinfection protocol.
3. Follow applicable manufacturer supplied instructions and warnings for any decontaminants, cleaning agents or other chemicals used to process ESTECH instruments.
4. Avoid prolonged instrument contact with saline solution, it can lead to pitting and stress corrosion cracking.
5. Do not allow blood or tissue to dry on instrument. Wash promptly after use. Long intervals between instrument use and processing for reuse (e.g., overnight or over the weekend) should be avoided to prevent contamination, corrosion and ineffective cleaning.
6. Avoid tap water. Water composition can have an adverse effect both on the treatment process and appearance of instruments and materials. Distilled water or deionized water for washing and sterilization are recommended.
7. Avoid using strongly acidic or alkaline detergents; use a neutral detergent. High alkaline solutions ($\text{pH} \geq 7$) can cause water marks and color changes in aluminum (AL6061, AL6060) instruments. High alkalines will increase the frequency of adjustment during preventive maintenance.
8. Avoid use of chlorines, mineral acids, solvents, phenol, iodine, peroxide and harsh abrasive agents by adhering to the following practices to prevent instrument rust and corrosion.
9. Avoid use of aldehyde organic compounds (i.e., glutaraldehyde, formaldehyde, etc.), which may cause blood soil to set on instrument.
10. Do not polish the surfaces of the product with visually coarse polishing powders, metal brush or metallic wools. They can limit the life of the instruments.
11. Do not overload instrument basket or tray.
12. All ESTECH surgical instruments with long, narrow cannulae, hinges and blind holes require particular attention during cleaning.
13. After washing, do not leave the product in a damp state.
14. No element of the process shall exceed the maximum temperature of 140°C.
15. Do not sterilize the product when any debris remains after washing.

LIMITATIONS ON REPROCESSING

1. Repeated processing has no significant influence on the length of life of the instruments.
2. End of life shall be determined by wear and damage from use.

EQUIPMENT NEEDED

1. Surgical scrub brush
2. Distilled water or deionized (DI) water
3. Lint-free cloth
4. Steam sterilizer
5. Ultrasonic cleaner (optional)
6. Instrument tray or basket
7. Enzymatic cleaner
8. Instrument Milk
9. *Geobacillus stearothermophilus* spores

FROM POINT OF USE

1. Clean new reusable instruments or as soon as possible after use. Wipe soil from instruments with a sponge or a towel moistened with water.
2. Soiled items should be stored in a manner to reduce the risk of exposure to pathogen and prevent damage to the instruments.
3. Place devices that were used in an instrument basket or tray to preclude blood and body fluids from drying on them. Cover instruments with a towel that has been moistened with water (not saline) or an enzymatic product specifically intended for stainless steel surgical instruments then place in an impervious bag. ESTECH recommends using ENZOL™ or a suitable equivalent. ENZOL is a brand of protein resolving cleaner to remove all traces of blood and debris.
4. Follow the manufacturer's instructions of the cleaning agent for the correct dilution/concentration, temperature and the correct water quality (e.g., pH, hardness) and correct exposure time.

DECONTAMINATION PREPARATION

1. Disassemble instruments only where intended without the use of tools. ESTECH shall provide, if applicable, specific instructions and tools for disassembling ESTECH medical devices.
2. Use flush port(s) to vigorously flush instrument with distilled water or deionized (DI) water or cleaning solution at least five (5) minutes at 22°C to 43°C (72°F to 110°F) to remove residual organic soil and lubricating agents without damaging the device.
3. As appropriate, loosen connections to allow brushes and cleaning solutions to flow in between parts.
4. Articulated instruments such as scissors, clamps, Universal Stabilizer Arms, forceps, etc. must be processed in the open position to minimize surface overlapping.

MANUAL CLEANING

1. The cleaning solution should be freshly prepared on a daily basis. When contaminate are high, it is recommended to prepare a fresh solution at even shorter intervals. Strictly follow the manufacturer's recommended instructions regarding temperature, concentration and proper use.
2. Immerse devices in the open position in an approved suitable enzymatic detergent instrument cleaner under recommended mixture concentration on the detergent's label for a minimum of 5 minutes to remove organic and inorganic materials. Note: Do not immerse an articulated arm (Rigid Table Mount Instrument Holder) in liquid.
3. Use a soft lint-free cloth or soft brush to manually scrub the device thoroughly (at least 1 minute) and parts thoroughly in intricate areas with enough friction to remove blood/bodily fluids, soil, tissue fragments, etc., from instrument.
4. Thoroughly clean all moving parts (e.g., beads, lumens, hinges, etc.) to prevent debris from interfering with movement.
5. Actuate and lightly tap instruments to release trapped air bubbles, check that the detergent diffuses to delicate areas of the device.
6. Rinse instruments thoroughly for a minimum of 1 minute with lukewarm distilled water or deionized (DI) water thoroughly to remove dislodged surgical debris and the detergent solution. Incomplete rinsing after cleaning can cause enzymatic cleaner to form a residue on the instrument.
7. Immediately after rinsing, gently drain off excess water. Blow dry using compressed air to blow off moisture. Dry instrument(s) thoroughly for minimum thirty (30) minutes).
8. As an alternate, wipe all the instruments with a clean, lint-free soft cloth repeatedly until completely dry.

ULTRASONIC WASHING

Note: It is recommended that the instruments and parts also be ultrasonically cleaned.

1. For ultrasonic treatment, instruments should be placed in open condition on proper perforated trays or in wire basket. Check that the device does not create "ultrasonic shadow areas" due to overlapping. Avoid any "wave shadows" or covering surfaces caused by wire baskets or perforated tray or bulky instruments.
2. Add *ENZOL* or an equivalent enzymatic cleaner *formulated for use by the manufacturer of the ultrasonic cleaner*. Strictly follow the manufacturer's instructions for concentration of the solution and proper use to avoid acidic or alkaline pH balances that can cause corrosion and breakage. Set the cycle times, suspension of the basket (e.g. the basket should not sit on the bottom of the ultrasonic cleaner), conditioning of the water, and other parameters according to the operating guide by the ultrasonic cleaner's manufacturer. Allow to sonicate for a minimum of 10 minutes.
3. Replace the detergent with a fresh supply periodically.
4. After ultrasonic treatment remove and rinse all instruments thoroughly with distilled or de-ionized water to avoid water spots. Rinse off all remnants of cleaning solution. Check that there are no loose parts. (e.g., screws, etc.) in the instrument cleaner.
5. Use compressed air to blow off moisture, wipe with a soft cloth, and dry completely for minimum thirty (30) minutes.

POST-CLEANING INSPECTION

Visually inspect the entire instrument (i.e., surface, beads, joints, air pockets, channels and lumen, etc.,) for cleanliness to assure dirt, soil and residue have been removed. If any tissue, blood, pus or soil is still present, repeat cleaning process. Check for corrosion, damaged surfaces and chipping. If product is damaged, it should not be used.

LUBRICATION

For optimal device performance and to maximize useful life, lubrication of the device using Steris® Hinge-Free® instrument lubricant (or equivalent) per manufacturer's instructions is recommended prior to each device sterilization cycle. Ensure the device is in an open position prior to lubrication.

STERILIZATION

After following the above cleaning recommendations, instruments may be sterilized using the following processing parameters. If applicable make sure instruments are not tightened during sterilization. To monitor the effectiveness of the sterilization process a biological indicator containing *Geobacillus stearothermophilus* spores may be used to verify a sterility assurance level (SAL) 10^{-6} kill of microorganisms.

Sterilizer Type	Temperature	Exposure Time
Gravity displacement, high-pressure steam sterilization, <i>unwrapped</i>	132°C (270°F)	3 minutes
Prevacuum (dynamic air-removal), <i>wrapped</i>	134°C (275°F)	3 minutes
Prevacuum (dynamic air-removal), <i>wrapped</i>	132°C (270°F)	4 minutes

POST-STERILIZATION INSPECTION

Prior to use inspect the entire instrument for proper function using the following verification test. If any wear is found, it is recommended the product be returned to ESTECH for replacement or repair. Do not use a device that is worn out, frayed, deteriorated, corroded, deformed, pitted, porous or otherwise damaged.

FUNCTIONAL CHECK

1. Perform a functional check for:
 - Damage and wear;
 - Cutting edges are free from nicks and present a continuous edge;
 - Accurate jaw and teeth alignment;
 - Smooth operation of moving parts thorough-out the intended range of motion;
 - Locking mechanisms are intact (close and open without difficulty); and
 - All components fit and assemble correctly with harmonizing components.

PREVENTIVE MAINTENANCE

1. To maintain optimum performance of the device, preventive maintenance once per year is recommended.
2. If the device is used frequently, or if a highly alkaline detergent with a pH 7 or higher is used, maintenance requirements may be more frequent.
3. When returning for maintenance or repair, please contact your local AtriCure sales representative or customer service for a RGA number. **Note:** Proper cleaning and decontamination is required prior to returning ESTECH instruments.

PACKAGING

Instruments must be packaged in a manner that will facilitate safe care and handling, e.g., non-biohazardous, protective container, tray with partitions; tips and edges protectors, etc. A certificate of decontamination must be included in the shipment.

WARRANTY

Instruments have a 1 year warranty on parts and workmanship. AtriCure warrants that the instruments are free from defects in both materials and workmanship. Suitability for use of the instruments for any surgical procedure shall be determined by the user. AtriCure shall not be liable for incidental or consequential damages of any kind. The above warranties are in lieu of all other warranties either expressed or implied including any warranty of any merchantability or fitness for use.

TECHNICAL SUPPORT


For more information about any ESTECH product or technical support, contact AtriCure customer service:

Manufacturer

Manufactured by:
AtriCure Incorporated
7555 Innovation Way
Mason, Ohio 45040 USA
Customer Service:
1-866-349-2342 (toll free)
1-513-755-4100 (phone)

European Authorized Representative

European Representative:
Köntges SPRL
Avenue Hellevelt 35
B-1180 Brussels Belgium
Tel: +32 (0) 2 375 51 63
FAX: +32 (0) 2 375 89 06
e-mail: herbert.kontges@skynet.be

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