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
**INSTRUCTIONS FOR USE** EN

**DESCRIPTION**

The AtriCure Dissector is a single patient use surgical instrument designed to dissect tissue during cardiothoracic surgical procedures. The Dissector’s battery-powered light source is used to navigate tissue for identification and isolation of anatomic structures.

**POWER SPECIFICATION**

- Voltage: 3.6 VDC
- Current: 35 mA
- Battery Capacity: 2100 mAh

**DISSECTOR ILLUSTRATION AND NOMENCLATURE (FIGURE 1)**

[1] Handle	[6] Rotation Knob
[2] Battery Door	[7] Light Source
[3] Battery Door Screw	[8] Dissecting Zone
[4] Pull-tab	[9] Hinge Point End
[5] Articulating End	[10] Suture Hole

**INDICATION FOR USE**

The AtriCure Dissector is intended to dissect cardiac tissue during cardiothoracic surgical procedures using AtriCure instrumentation.

**CONTRAINDICATIONS**

The Dissector is not a sealed device and is not intended for insertion through ports which maintain insufflation.

**⚠️ WARNINGS ⚠️**

To avoid the risk of patient infection, inspect the product packaging prior to opening to ensure that the sterility barrier is not breached. If the sterility barrier is breached, do not use the Dissector.

Do not immerse any part of the Dissector in liquids as this may damage the device.

Do not touch the Dissector tip against metal staples or clips. This may cause damage to the light source lens.

The use of other cables and accessories may affect the Dissector’s EMC performance.

DO not use excessive force when articulating the Dissector. Using excessive force when articulating may damage tissue or the Dissector hinge point.

Avoid contacting the Dissector with electrodes of any electro-surgical device. This may damage the Dissector, the electro-surgical device, or tissue.

The Dissector should not be used adjacent or stacked with other equipment.

Battery Warning. The Dissector contains a Lithium disposable battery. Do not recharge, disassemble, heat above 100°C, incinerate, or expose the battery directly to water.

**⚠️ CAUTIONS**

- Carefully read all instructions for the AtriCure Dissector prior to using the device. Failure to properly follow instructions may lead to injury and/or may result in improper functioning of the device.
- Medical electrical equipment needs special precautions regarding EMC and needs to be installed according to EMC information.
- Use of the Dissector should be limited to properly trained and qualified medical personnel.
- The Dissector can radiate radio frequency energy and, if not installed, used, and serviced in accordance with electromagnetic compatibility information provided in the instructions, may cause harmful interference to other devices in the vicinity.
- Portable and mobile RF communications equipment can also affect Dissector performance and care should be taken to minimize such interference. However, there is no guarantee that interference will not occur in a particular installation.
- Variations in specific procedures may occur due to individual physician techniques and patient anatomy.
- To avoid damage to the device, do not drop or toss the Dissector. If the Dissector is dropped, do not use. Replace with a new Dissector.
- To avoid damage to the Dissector, use caution when handling and interfacing with other devices.
- Do not re-sterilize or reuse the Dissector. Single Patient use only. Reuse can cause patient injury and/or the communication of infectious disease(s) from one patient to another.
- Do not look directly into the light source when illuminated.
- During a surgical procedure, ensure that the Dissector hinge point remains visible at all times. The hinge point should always be visible for the frame of reference of the tip location.
- The dissector includes an LED light source intended to indicate the device position and orientation, not to support visualization of structures. Appearance shifts occur during the use of LED light sources due to the difference in color, temperature, and CRI characteristics of the LED light from those of normal white light sources.
- The Dissector shall operate throughout ambient temperature, humidity and atmospheric pressure ranges consistent with a surgical operating room.
- The Dissector is intended to be used in the surgical operating room.
- The Dissector will be used near HF SURGICAL EQUIPMENT.
- The Dissector is not intended for use in the RF shielded room of an ME SYSTEM for magnetic resonance imaging, where the intensity of EM DISTURBANCES is high.
- Be aware of that other electrical devices may produce electromagnetic interference which could impact the performance of this device. In such instances remove the source of the interference or increase the separation between the source of interference and this device. The Dissector has no essential performance that can be affected by EM DISTURBANCES.

**INSTRUCTIONS FOR USE**

- Using sterile technique, remove the Dissector from its packaging.
- Remove the pull-tab from the Dissector handle to activate the light source. Discard the pull-tab and ensure that the light source is illuminated. Do not attempt to reinsert pull-tab.
- The Dissector will automatically stay illuminated for a minimum of 4 hours of useful life. After 4 hours, the Dissector may remain illuminated until all power is discharged or the battery is removed.
- To verify the range of motion, rotate the articulation knob before inserting into the operating field.
  - The handle provides the user with an understanding of device orientation. The top of the device handle includes a flat surface which indicates up and orients the plane of dissection. The bottom of the device handle includes a curved surface with finger grips.
  - Rotating the knob counterclockwise articulates the distal end away from the user.
  - Rotating the knob clockwise articulates the distal end toward the user. The Dissector tip should move freely in either direction.
- Before each use, manually check the Dissector to ensure that there are no sharp edges. If there are any sharp edges, do not use the device. Replace with a new Dissector.
- Using the rotation knob, adjust the tip so it is positioned away from the user.
- Insert the Dissector into the targeted dissection plane.

**NOTE:** A tissue incision of 12 mm (0.47 inches) minimum is recommended for insertion of the Dissector.

- During use, ensure that the articulation hinge remains exposed. The articulation hinge should be used to guide the depth and orientation of the dissecting tip.
  - Be sure to maintain support of the Dissector while in the surgical field.
- Using the light source to obtain placement feedback, articulate the rotation knob as needed to facilitate blunt dissection of the desired tissue.
  - A soft tipped tool may be used to assist in tissue dissection.
  - The surgeon must choose the appropriate guide.
- Inspect the surgical area to ensure that dissection is achieved.
- Make appropriate adjustments to the rotation knob to assist in removal of the Dissector.
- Discard the Dissector after use. Follow the local governing ordinances and recycling plans regarding disposal or recycling of device components. Remove the Lithium battery prior to discarding the Dissector. The battery may be discarded by removing the battery door screw, opening the door, and disconnecting the battery. Dispose of the battery in an appropriate manner.

**Table 1: Guidance and Manufacturer’s Declaration – Electromagnetic Emissions**

The AtriCure Dissector is intended for use in the electromagnetic environment specified below. The customer or the user of the Dissector should assure that it is used in such an environment.		
Emissions Test	Compliance	Electromagnetic Environment – Guidance
RF Emissions CISPR 11	Group 1	The Dissector uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class B	The Dissector is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic Emissions IEC 61000-3-2	Not Applicable	
Voltage Fluctuations/Flicker Emissions IEC 61000-3-3	Not Applicable	

**Table 2: Guidance and Manufacturer’s Declaration – Electromagnetic Immunity**

The AtriCure Dissector is intended for use in the electromagnetic environment specified below. The customer or the user of the Dissector should assure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	Not applicable	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical Fast Transient / Burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV line(s) to line(s) ±2 kV line(s) to earth	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Voltage Dips, Short Interruptions and Voltage Variations on Power Supply Input Lines IEC 61000-4-11	Voltage Dips: 0% U <sub>i</sub> ; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° phase angles 0% U <sub>i</sub> ; 1 cycle and 70% U <sub>i</sub> 25/30 cycles Single phase: at 0° Voltage interruptions: 0% U <sub>i</sub> ; 250/300 cycle	Not applicable	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Dissector requires continued operation during power mains interruptions, it is recommended that the Dissector be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) Magnetic Field IEC 61000-4-8	30 A/m	Not applicable	Power frequency magnetic fields should be at levels characteristic of a typical commercial or hospital environment.

NOTE: U<sub>i</sub> is the a.c. mains voltage prior to application of the test level.

**HOW SUPPLIED**

The Dissector is supplied STERILE and NON-PYROGENIC in unopened, undamaged package. For single use only. Do not re-sterilize. Do Not Re-Use.

**RETURN OF USED PRODUCT**

If for any reason this product must be returned to AtriCure, Inc., a return goods authorization (RGA) number is required from AtriCure, Inc., prior to shipping.

If the product has been in contact with blood or body fluids, it must be thoroughly cleaned and disinfected before packing. It should be shipped in either the original carton or an equivalent carton, to prevent damage during shipment; and it should be properly labeled with an RGA number and an indication of the biohazardous nature of the contents of shipment.

Instructions for cleaning and materials, including appropriate shipping containers, proper labeling, and an RGA number may be obtained from AtriCure, Inc.

**⚠️ CAUTION:** It is the responsibility of the health care institution to adequately prepare and identify the product for shipment.

**DISCLAIMER STATEMENTS**

Users assume responsibility for approving the acceptable condition of this product before it is used, and for ensuring that the product is only used in the manner described in these instructions for use, including, but not limited to, ensuring that the product is not re-used.

Under no circumstances will AtriCure, Inc. be responsible for any incidental, special or consequential loss, damage, or expense, which is the result of the deliberate misuse or re-use of this product, including any loss, damage, or expense which is related to personal injury or damage to property.

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