## Sensing Cable Kit: CSK-2010 Instructions for Use

 AtriCure Incorporated  
7555 Innovation Way  
Mason, Ohio 45040 USA  
Customer Service:  
1-866-349-2342 (toll free)  
1-513-755-4100 (phone)

 AtriCure Europe B.V.  
De entree 260  
1101 EE Amsterdam  
NL  
+31 20 7005560  
ear@atricure.com

For Use only with the nContact Coagulation System

### INDICATIONS:

For use only with nContact Generator.  
The CSK-2010 Sensing Cable is provided non-sterile and should not be sterilized. See cleaning instructions section of this document.

### CONTENTS:

**CSK-2010 Sensing Cable Kit contains:**

1 ea CS-2010 Sensing Cable                      1 ea Instructions for Use (IFU) Brochure

### INSTRUCTIONS:

#### CSK-2010

1. Remove cable from package.
2. Before every use inspect cable for damage. Pay particular attention to the insulation of the cable. If any part of the cable is damaged, discard without using. See next page for schematic.
3. Connect GREY end of Sensing Cable to the GREY port “Diagnostic” on nContact Generator until a click is felt.
4. Connect BLACK end of CSK-2000 Cable to the BLACK receptacle of the generator “Coagulation” until a click is felt.
5. Connect 2mm shrouded pins to EP sensing equipment.
6. Follow instructions for use of generator and coagulation device.
7. To disconnect Sensing Cable from the generator, grasp by connector cover and pull back. Do NOT pull on wire.
8. To disconnect 2mm shrouded pins, grasp the over-molded portion of the 2mm pin where the pin number is located and pull. Do NOT pull on the wire.
9. Excess pulling and flexing of the cable may damage the cable rendering it non-functional.
10. Store cable in safe dry area until next use.

### CLEANING INSTRUCTIONS:

 CAUTION: Ensure Isopropyl Alcohol (IPA) is completely dry before operating the unit.

 CAUTION: Avoid caustic or abrasive cleaners

1. Disconnect the cable from the unit before cleaning.
2. If the cable is contaminated with blood and other body fluids, it shall be cleaned before the contamination can dry (within two hours of contamination).
3. The outer surfaces of the cable shall be cleaned with 70% - 90% Isopropyl Alcohol (IPA) wipes for a minimum of two minutes.

4. Pay attention to all area where fluids or soil may gather, such as under/around the connector or any tight crevices/grooves.
5. Dry the cable with a dry, white lint-free cloth.
6. Conduct a final confirmation of the cleaning process by visually inspecting the white cloth for remaining soil.
7. If soil remains on the white cloth, repeat steps 3 through 6.
8. After cleaning is complete, turn the unit on to perform Power On Self-Test (POST). If any errors are received, contact AtriCure to begin return process.

### DISPOSAL:

Follow local governing ordinances and recycling plans regarding disposal or recycling of device components.

### PRECAUTIONS:


1. Inappropriate handling of the connecting cable, including sterilization and immersing the electrical connectors, can result in degradation of system performance including the inability to initiate or complete and ablation therapy.

### WARNINGS:

- For information concerning all warnings, precautions, and troubleshooting see nContact Coagulation System Radiofrequency (RF) Generator Unit Model CS-3000 Operators Manual.
- Ensure that the EP Sensing (EKG) equipment comply with IEC 60601-2-25 for protection from high frequency surgical interference.
- If necessary, clean the cable with 70%- 90% Isopropyl Alcohol (IPA) wipes to ensure that the cable does not cause an infection.

 **CAUTION:** Cables to surgical electrodes should be positioned to prevent contact with patient or other leads.

### EXPLANATION OF SYMBOLS FROM PACKAGE LABEL:

	Caution		Lot Number
	Non-sterile reusable device		Model Number
	CE Mark and Identification number of Notified Body		Date of Manufacture
	Follow instructions for use		Manufacturer
	Temperature Storage Limit		Humidity Storage Limit
	Federal (US) law restricts this device to sale by or on the order of a physician.		Waste Electrical Electronic Equipment
	Authorized Representative in the European Community		

### NOMENCLATURE:


RF=Radiofrequency                      LBL=Label

### WARRANTY INFORMATION:

AtriCure warrants that reasonable care has been used in the design and manufacture of this instrument. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise,

including, but not limited to, any implied warranties of merchantability or fitness for a particular use. Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond AtriCure control directly affect the instrument and the result obtained from its use. AtriCure's obligation under this warranty is limited to the repair or replacement of this instrument and AtriCure shall not be liable for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of this instrument. AtriCure neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. AtriCure assumes no liability with respect to instruments reused, reprocessed or re-sterilized and makes no warranties expressed or implied, including but not limited to merchantability or fitness for intended use, with respect to such instrument.

Under no circumstances will AtriCure, Inc. be responsible for any incidental, special or consequential loss, damage, or expense, which is the result of the deliberate misuse or re-use of this product, including any loss, damage, or expense which is related to personal injury or damage to property.

	English	Nederlands	Deutsch	Français	Italiano	Español	Português
<b>RF</b>	RF Coil Electrode	RF-spoelelektrode	RF Spulenelektrode	Électrode de la bobine RF	Elettrodo bobina RF	Electrodo de bobina de RF	Elétrodo de Bobina de RF
<b>D1</b>	D1 Sensing Electrode	D1-sensingelektrode	D1 Sensorelektrode	Électrodes de détection D1	Elettrodo di rilevamento D1	Electrodos de detección D1	Elétrodo de Medição D1
<b>D2</b>	D2 Sensing Electrode	D2-sensingelektrode	D2 Sensorelektrode	Électrodes de détection D2	Elettrodo di rilevamento D2	Electrodos de detección D2	Elétrodo de Medição D2
<b>P1</b>	P1 Sensing Electrode	P1-sensingelektrode	P1 Sensorelektrode	Électrodes de détection P1	Elettrodo di rilevamento P1	Electrodos de detección P1	Elétrodo de Medição P1
<b>P2</b>	P2 Sensing Electrode	P2-sensingelektrode	P2 Sensorelektrode	Électrodes de détection P2	Elettrodo di rilevamento P2	Electrodos de detección P2	Elétrodo de Medição P2
<b>#1=D1</b>	Shrouded Pin #1, D1 Sensing Electrode	Omhulde pin nr.1, D1-sensingelektrode	Ummantelter Stift Nr. 1, D1 Sensorelektrode	Broche protégée n°1, électrode de détection D1	Pin schermato #1, elettrodo di rilevamento D1	Perno cubierto n.º 1, electrodo de detección D1	Pino Blindado N°1, Elétrodo do Sensor D1
<b>#2=D2</b>	Shrouded Pin #2, D2 Sensing Electrode	Omhulde pin nr.2, D2-sensingelektrode	Ummantelter Stift Nr. 2, D2 Sensorelektrode	Broche protégée n°2, électrode de détection D2	Pin schermato #2, elettrodo di rilevamento D2	Perno cubierto n.º 2, electrodo de detección D2	Pino Blindado N°2, Elétrodo do Sensor D2
<b>#3=P1</b>	Shrouded Pin #3, P1 Sensing Electrode	Omhulde pin nr.3, P1-sensingelektrode	Ummantelter Stift Nr. 3, P1 Sensorelektrode	Broche protégée n°3, électrode de détection P1	Pin schermato #3, elettrodo di rilevamento P1	Perno cubierto n.º 3, electrodo de detección P1	Pino Blindado N°3, Elétrodo do Sensor P1
<b>#4=P2</b>	Shrouded Pin #4, P2 Sensing Electrode	Omhulde pin nr.4, P2-sensingelektrode	Ummantelter Stift Nr. 4, P2 Sensorelektrode	Broche protégée n°4, électrode de détection P2	Pin schermato #4, elettrodo di rilevamento P2	Perno cubierto n.º 4, electrodo de detección P2	Pino Blindado N°4, Elétrodo do Sensor P2
<b>A</b>	TO SENSING or EKG EQUIPMENT	NAAR SENSING-of ECG-APPARATUUR	ZUR SENSOR- oder EKG-VORRICHTUNG	VERS L'ÉQUIPEMENT DE DÉTECTION ou EKG	ALL'APPARECCHIATURA DI RILEVAMENTO O EKG	A EQUIPOS ELECTRO-FISIOLÓGICOS O DE DETECCIÓN	PARA EQUIPAMENTO SENSOR ou ECG

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