

**U.S. Indications: For Adult Patients** AtriCure's cryoICE® cryoSPHERE® cryoablation probes are sterile, single use devices intended for use performed by freezing target tissues, creating an inflammatory response (cryonecrosis) for blocking pain by temporarily ablating peripheral nerves. **For Adolescent Patients** The cryoICE cryoSPHERE cryoablation probes are intended for use to temporarily block pain by ablating intercostal nerves under direct visualization<sup>1</sup> in adolescent patients of at least 12 years of age.

<sup>1</sup>Direct visualization, in this context, requires that the surgeon is able to see the targeted tissue for cryoablation directly or with assistance from a camera, endoscope or other similar optical technology.

Please review the Instructions for Use for a complete listing of contraindications, warnings, precautions and potential adverse events prior to using these devices.

**Rx Only.**

**ATRICURE, INC.**  
7555 Innovation Way  
Mason, Ohio 45040 USA  
Customer Service: +1 (866) 349-2342  
CustomerService@AtriCure.com  
[www.AtriCure.com](http://www.AtriCure.com)

PM-US-0058E-0425-G

# *cryo*SPHERE®

A non-pharmacological/  
non-opioid solution for  
post-operative pain.



## *cryo*SPHERE®

- Enhanced Flexibility
- Uniform Application
- Engineered for Intercostal Nerve Cryoablation


### Indications

**For Adult Patients:** AtriCure's cryoICE® cryoSPHERE® cryoablation probes are sterile, single use devices intended for use in blocking pain by temporarily ablating peripheral nerves.

**AtriCure**

**AtriCure**

**A non-pharmacological/non-opioid solution that can provide effective post-operative pain management for several months.**



**Features**

**8mm Ball Tip**      Designed for minimally invasive applications and placement on intercostal nerve

**Bendable Distal Shaft**      Bendable distal shaft up to 180 degrees


**Braided Insulation**      Technology that is designed to maintain a temperature warmer than -10°C in tissue

**Reinforced Proximal Shaft**      Provides rigidity at the level of the incision points

**Ergonomic Handle Design**      Provides ease of use and comfort during cryoablation

**Active Defrost**      Allows probe to be easily removed from cryoablation site

**Malleable Cryoablation Probe**



**cryoICE BOX® - V6**

ACM1 (US)



**cryoICE Cart & Upgrade Kits**

ASCS (cart with double tank holder, US only)

ASCC (single tank holder upgrade)

ASCC2 (double tank holder upgrade)

20 lbs N<sub>2</sub>O tanks not provided by AtriCure



**Foot Switch (optional)**

For use with ACM1/ACM2  
FTSWITCH

