

AtriClip PRO2 Device

AtriClip PRO Devices can be used in totally thoracoscopic procedures in conjunction with other heart heart surgery.



Features

- Hoopless End Effector
- Small Diameter End Effector (12 mm)
- Active Articulation Levers

| LAA Exclusion System | | |
|----------------------|---------------------|----------------------|
| Size | AtriClip PRO Device | AtriClip PRO2 Device |
| 35 mm | PRO135 | PRO235 |
| 40 mm | PRO140 | PRO240 |
| 45 mm | PRO145 | PRO245 |
| 50 mm | PRO150 | PRO250 |

U.S. Indications: The AtriClip LAA Exclusion System is indicated for the occlusion of the heart's left atrial appendage, performed under direct visualization and in conjunction with other cardiac surgical procedures.

Direct visualization, in this context, requires that the surgeon is able to see the heart directly, with or without assistance from a camera, endoscope, etc., or other appropriate viewing technologies.

Please review the Instructions for Use for a complete listing of contraindications, warnings, precautions and potential adverse events prior to using these devices.

Rx Only.

ATRICURE, INC.
7555 Innovation Way
Mason, Ohio 45040 USA
Customer Service: +1 (866) 349-2342
CustomerService@AtriCure.com
www.AtriCure.com

AtriClip PRO2[®] Device

A Simplified Approach to LAA Exclusion




AtriClip PRO2 Device

AtriClip PRO[®]
Devices can now
be used in totally
thoracoscopic
procedures.



Greater than 30°
Pitch and Yaw

Rigid Shaft



12 mm
Diameter
End
Effector

